Invia[®] Motion[™]

NEGATIVE PRESSURE WOUND THERAPY SYSTEM

- EN Clinician instructions for use
- DE Gebrauchsanweisung für medizinisches Personal
- FR Mode d'emploi à l'usage du personnel médical
- ES Instrucciones de uso para profesionales sanitarios
- NL Gebruiksaanwijzing voor artsen
- SV Bruksanvisning för vårdpersonal
- IT Istruzioni per l'uso per il medico
- FI Käyttöohjeet lääkärille
- PT Instruções de utilização para o médico
- DA Brugsanvisning til sundhedspersonalet
- NO Bruksanvisning for helsepersonell
- SK Návod na používanie pre lekárov

Table of contents

| | |
|--|----|
| Introduction | 4 |
| Overview | 4 |
| Intended use | 6 |
| Indications for use | 6 |
| Contraindications..... | 6 |
| Warnings and safety instructions | 7 |
| Warnings | 7 |
| Cautions | 8 |
| Physician's orders..... | 9 |
| Wound assessment..... | 9 |
| Dressing technique..... | 10 |
| Initial delivery | 11 |
| Invia Motion NPWT pump accessories and associated products | 11 |
| Invia Motion disposables | 12 |
| Wound dressings..... | 12 |
| Charge battery for the first time..... | 13 |
| Disconnect from power source..... | 13 |
| Positioning..... | 13 |
| Preparation for use | 14 |
| Operating instructions | 17 |
| Power supply and battery management | 19 |
| Notifications and troubleshooting..... | 21 |
| Overview notifications..... | 24 |
| Cleaning and disinfection..... | 25 |
| Disposal | 26 |
| Warranty and maintenance..... | 27 |
| Signs and symbols | 28 |
| Technical specifications | 30 |
| Electromagnetic compatibility (EMC)..... | 31 |

Introduction

The Invia Motion Negative Pressure Wound Therapy system, subsequently named "Invia Motion NPWT system", is approved exclusively for the use described in these Instructions for use. Please read the information contained within these Instructions for use before operation of the Invia Motion NPWT system. These instructions for use must be kept with the device.

These instructions for use are a general guide for the use of the product. Medical situations must be addressed by a physician.

Indications, Contraindications, Warnings, Cautions and other information to support the safe use of the Invia Motion NPWT system are included in the respective chapters in these instructions for use.

The Invia Motion NPWT pump is a suction pump for Negative Pressure Wound Therapy that provides treatment status through a display and acoustic signals. The Invia Motion NPWT pump is a single patient use pump which provides continuous or intermittent operation and multiple negative pressure selection options. The Invia Motion NPWT pump is portable and can be operated using a rechargeable battery. Acoustic and optical signals are triggered for variances from the set values as well as for faults. The Invia Motion NPWT system is intended for use in all healthcare settings, including home care, by healthcare personnel and caregivers trained on the instructions for use. The patient or caregiver must have the visual and hearing acuity necessary to appropriately respond to notifications with the sensory and cognitive ability to understand directions from a healthcare professional pertaining to the proper use of the pump.


Invia Motion Negative Pressure Wound Therapy system is available in two versions:

- 1) **REF** 087.4022: Invia Motion – 15 Days Negative Pressure Wound Therapy system
- 2) **REF** 087.4024: Invia Motion – 60 Days Negative Pressure Wound Therapy system

Overview

The Invia Motion NPWT system and the Invia Motion carrying case are for a single patient's entire therapy.

Buttons to operate the Invia Motion NPWT pump

Backlight

The backlight will be activated by pressing any button. The backlight will turn off after 10 seconds if no further interaction results.

Symbols on display

| Notifications | |
|---------------|---|
| | Leakage |
| | Blockage |
| | Canister full |
| | Battery low (symbol is blinking) |
| | Battery fully discharged (symbol is blinking) |
| | Battery missing (symbol is blinking) |
| | Internal fault |
| | Pump time expired (symbol is blinking) (Example for Invia Motion – 15 Days) |
| Information | |
| | Shows runtime already used in days (📅) and hours (🕒) (Example for 6 days and 23 hours) |
| | Pump is running / On-state in intermittent mode (symbol is blinking) |
| | Off-state in intermittent mode (symbol is blinking) |
| | Constant mode |
| | Intermittent mode |
| | Negative pressure setting in mmHg |
| | Battery fully charged |
| | Battery charge level |
| | Battery charging (bars moving) |

Intended use

The Invia Motion NPWT system is intended to be used by healthcare professionals or adequately trained lay users.

Healthcare professionals are responsible to train lay users according to the patient instructions for use and explain all related safety information.

Intended patient population

The Invia Motion NPWT system is intended to be used on patients only exhibiting conditions as described in the indications for use. The device has not been studied in pediatric patients.

Intended environment

The Invia Motion NPWT system is intended for use in acute, extended and home care settings.

Indications for use

The Invia Motion Negative Pressure Wound Therapy (NPWT) system is indicated for patients who would benefit from a suction device (NPWT) as when used on open wounds it creates an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. When used on closed surgical incisions, the Invia Motion NPWT system is also intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudate via the application of Negative Pressure Wound Therapy.

The Invia Motion NPWT system is appropriate for the following indications:

- Acute or subacute wounds
- Chronic wounds
- Dehisced wounds
- Pressure ulcers
- Diabetic/Neuropathic ulcers
- Venous insufficiency ulcers
- Traumatic wounds
- Partial thickness burns
- Flaps and grafts
- Closed surgical incisions

Contraindications

The Invia Motion NPWT system is contraindicated in the presence of:

- Necrotic tissue with eschar present
- Untreated osteomyelitis
- Non-enteric and unexplored fistulas
- Malignancy in the wound
- Exposed vasculature
- Exposed nerves
- Exposed anastomotic site of blood vessels or bypasses
- Exposed organs

Warnings and safety instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety related tip

Indicating useful information about the safe use of the device.

The Invia Motion Negative Pressure Wound Therapy system is intended for use as described in these instructions for use.

Medela is only responsible for the effect on BASIC SAFETY, reliability and performance of the Invia Motion NPWT system if it is used in accordance with the instructions for use.

**Please read and observe these warnings and safety instructions before operation.
These instructions for use must be kept with the device.**

Please note that these instructions for use are a general guide for the use of the product. Medical situations must be addressed by a physician.

Warnings

WARNINGS

- Do not modify this equipment without authorization from the manufacturer.
- Consult the indications for use, contraindications, warnings and cautions before using the Invia Motion NPWT pump as a negative pressure source with the Invia foam and gauze dressings. Failure to read and follow all instructions in this manual prior to use may result in serious or fatal injury of the patient.
- Invia wound therapy Instructions advise 24 hours therapy per day without interruption. If therapy is discontinued for more than 2 hours, the dressing should be replaced and therapy restarted by a healthcare professional.
- Unplug the pump before cleaning.
- Do not use accessories other than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased electromagnetic emissions or decreased electromagnetic immunity of the Invia Motion NPWT pump.
- The Invia Motion NPWT system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Motion NPWT system should be observed to verify normal operation in the configuration in which it will be used.
- The Invia Motion NPWT pump is intended for a single patient's entire therapy. Reuse on another patient might cause cross contamination.
- Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately stop use of the pump, apply pressure and seek immediate emergency medical attention.
- Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger to the patient.
- When using the Invia Motion NPWT system in close proximity to blood vessels, organs, nerves, tendons or ligaments a protective barrier, such as a non-adherent contact layer, must be used.
- Serious or fatal injury can result from bone fragments or sharp edges (e.g. staples or hardware) that could puncture protective barriers, vessels or organs.

- Patients with spinal cord injuries (stimulation of sympathetic nervous system) who experience Autonomic Hyperreflexia during therapy should have NPWT discontinued.
- In the event that defibrillation is required, disconnect the pump from the wound dressing before the patient is defibrillated.
- Special care is advised for dressing placement and removal in order to avoid situations such as unintentional gauze or foam retention.
- This device has not been studied in pediatric patients.
- Do not bring the device into the following areas: Positron Emission Tomography (PET) or Hyperbaric Oxygen Chamber (HBO).
- The Invia Motion NPWT system is Magnetic Resonance (MR) Unsafe. Do not take the Invia Motion NPWT system into the MR environment.
- HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Invia Motion pump.
- Explosion hazard – the Invia Motion NPWT pump is not for use in potentially explosive environments including oxygen enriched environments and in areas of flammable anesthetics.
- Portable Radio Frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 0.3 m to any part of the Invia Motion NPWT system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- The Invia Motion NPWT system is intended to be used in conjunction with the Invia foam and gauze dressings.
- Invia Foam Dressing Kit with FitPad, X-Large should not be used in conjunction with the Invia Motion NPWT pump.

Cautions

CAUTIONS

- Incorrect use may cause pain and injury to the patient. Excessive negative pressure or an infection of the wound may cause pain to the patient. In either case, the dressing must be changed and the wound inspected in detail.
- The patient should be monitored regularly according to physician instructions and facility guidelines to monitor patient comfort, therapy compliance and signs of wound infection.
- Do not use an Invia Motion canister/tubing set if the sterile packaging is damaged.
- Wireless communication equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, and walkie-talkies can interfere with the Invia Motion pump and should be kept at a minimum distance of 30 cm away from the Invia Motion pump.
- The Invia Motion NPWT pump must remain in an upright (carrying case) or lying (display side up) position during use.
- The total wound volume connected to the Invia Motion NPWT pump shall be 1l or less.

The following statements describe medical conditions that may require special care for the safe and effective use of the Invia Motion NPWT system.

- Patients at high risk for bleeding and hemorrhage, including patients experiencing active bleeding or difficult wound hemostasis.
- Patients taking anticoagulants or platelet aggregation inhibitors.
- Patients with a history of vascular anastomosis or friable, irradiated, sutured or infected blood vessels.
- Patients with spinal cord injury (sympathetic nervous system stimulation): in the event a patient experiences autonomic dysreflexia (sudden changes in blood pressure or heart rate in response to stimulation of the sympathetic nervous system), discontinue negative pressure wound therapy to help minimize sensory stimulation and seek immediate medical assistance.

- Patients with infected wound or osteomyelitis.
- Wounds that involve an enteric fistula.
- To minimize the risk of bradycardia, NPWT must not be placed in the proximity of the vagus nerve.
- Avoid circumferential dressing applications.
- The Negative Pressure Wound Therapy must be used 24 hours per day without interruption. If the pump is stopped for more than 2 hours, the dressing must be changed and therapy restarted.
- Consider the patient's size and weight when prescribing this device.
- Consider mode of therapy – intermittent versus continuous.
- For maximum benefit on closed surgical incision, the Invia negative pressure therapy should be applied immediately post-surgery to clean surgically closed wounds. It is to be continuously applied for a minimum of two days up to a maximum of seven days, with regular dressing changes every 48 to 72 hours. All dressing changes should be applied under direct medical supervision. The Invia therapy system will not be effective in addressing complications associated with the following:
 - Ischemia to the incision or the incision area
 - Untreated or inadequately treated infection
 - Inadequate hemostasis of the incision
 - Cellulitis of the incision area

Physician's orders

A medical professional must assess each wound when applying this therapy and determining negative pressure intensity according to the wound characteristics. Therapy or setting changes may only be done under the order of a physician.

Wound assessment

CAUTION

Patient monitoring: The patient should be monitored regularly according to the physician's instructions and facility guidelines to check for patient comfort, therapy compliance and signs of infection.

WARNING

Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger to the patient.

Observe the wound/periwound tissue and exudate for signs of infection or other complications. The most common signs of infection include redness, tenderness, fever, swelling, itching, increased warmth in the wound area, strong odor or purulent discharge. Additional symptoms include nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucous membranes, disorientation, high fever (>38.8° C), refractory hypotension, orthostatic hypotension, or erythroedema (a sunburn-like rash). More serious complications of infection include pain, discomfort, fever, gangrene, toxic or septic shock. If more serious complications of infection occur, discontinue therapy and consult a healthcare professional immediately.

Dressing technique

 Wound dressings to be applied and changed by healthcare professionals only.

Consult the appropriate Invia dressing instructions for use for information regarding dressing applications included in every Invia dressing carton. Perform a thorough wound cleansing per physician orders prior to dressing applications.

 For closed surgical incision the recommended pressure level is -125 mmHg at constant mode.

Routine dressing changes should occur every 48–72 hours. Dressing changes for infected wounds should be considered more frequently.

Initial delivery

Upon delivery, check the package containing the Invia Motion NPWT system for completeness and general condition.

1 Invia Motion NPWT pump

1 Invia Motion carrying case

1 Invia Motion power supply

1 Invia Motion Clinician Instructions for Use
1 Invia Motion Patient Instructions for Use
1 Invia Motion Quick Card

☞ When first turned on, the Invia Motion NPWT pump needs to be attached to a wall outlet in order to fully charge the battery.

Invia Motion NPWT pump accessories and associated products

The Invia Motion NPWT system was verified in combination with the accessories listed below. For correct and safe operation use the Invia Motion NPWT system with these accessories only:

Y-connector
REF087.0023

Invia Motion NPWT pump

Power supply
REF101033647

Carrying case
REF087.4003

Canister/Tubing Set
REF087.4016

Wound dressings

 Wound dressings to be applied and changed by healthcare professionals only.

The Invia Motion NPWT system is intended to be used in conjunction with the Invia dressings only. For specific dressing indications, contraindications, warnings and cautions, consult the appropriate Invia dressing instructions for use.

WARNING

Invia foam dressing kit with FitPad, X-Large should not be used in conjunction with the Invia Motion NPWT pump.

CAUTION

The total wound volume connected to the Invia Motion NPWT pump shall be 1l or less.

Invia foam and gauze dressing kits provide their own instructions for use.


Invia Motion disposables

Invia Motion canister/tubing set 0.15l (150 ml), with solidifier (sterile, single use)

Capacity: 0.15l (150 ml)
Material: Polypropylene, Silicon, ABS
Sterilization Method: Ethylene Oxide

Intended use


The Medela Invia Motion canister/tubing set 0.15l (150 ml) is an accessory of the Medela Invia Motion NPWT system and it is intended to collect removed wound exudate and infectious material.

Charge battery for the first time

Fully charge battery before first use. Charging time can be up to 4 hours. Select the plug. Attach the plug to the power supply.

Plug the Invia Motion power supply into a wall outlet and open the small protective cover of the pump for the Direct Current (DC) socket and attach the power supply. Moving bars on the battery display will indicate charging is occurring. [].

1. Battery is fully discharged

2. Battery is being charged (moving bars)

3. Battery is fully charged

Disconnect from power source

 Check the battery charge level periodically. Make sure the battery compartment on the back side of the pump is closed. Do not remove battery from device as long as therapy is ongoing.

Disconnecting the power source

When the battery is fully charged (charging time can be up to 4 hours) the following non-blinking symbol will appear on the display . Charging of the battery can be slowed when exposed to high temperatures. To disconnect the Invia Motion NPWT pump from the electrical source, hold the end of the plug and remove it from the pump.

Close the DC socket by pressing the small protective cover over the opening.

Positioning

CAUTION

The Invia Motion NPWT pump must remain in an upright (carrying case) or lying (display side up) position during use.

The Invia Motion NPWT pump is protected against the penetration from solid/fluid substances by a hydrophilic bacteria filter.

If the canister is full or the Invia Motion NPWT pump is placed in an inappropriate position, fluids will flow into the safety chamber and clog up the filter. In this situation, the pressure is interrupted and canister full notification will appear (acoustic signal and symbol on display). The Invia Motion canister/tubing set must be replaced.

Preparation for use

 To be used by healthcare professionals or adequately trained caregivers only.

Checklist

Prior to using the Invia Motion NPWT system, it is important to check for the following:

- Completeness and general condition of the packaging
- Damage to the power cord and plug
- Damage to the pump
- Proper functioning of the device
- Battery status

Connect Invia Motion canister/tubing set 0.15 l (150 ml)

 Make sure the correct canister is used: **REF** 087.4016.

1. Open sterile packaging.

CAUTION

Do not use the Invia Motion canister/tubing set if the sterile packaging is damaged or has expired.


2. Insert the Invia Motion canister/tubing set into the Invia Motion NPWT pump connecting port (1). Rotate the canister/tubing set to the right to secure (2).

3. Make sure the canister is aligned with the label to the front.

4. Connect the dressing tubing to the pump tubing.

Make sure both parts are aligned correctly as shown in the picture below. Push together until you hear a click.


For dressing applications, please refer to the instructions for use provided with the Invia dressing kits.

Switch Invia Motion NPWT pump on

1. Press [] to switch the Invia Motion NPWT pump on.

2. When turned on for the first time, the Invia Motion NPWT pump will start up in the default factory setting (constant mode with pressure level of -125mmHg).

If the Invia Motion NPWT pump is turned off, the setting will be as previously set when the pump is turned back on.


Run-time, timer and safety time

REF 087.4022 Invia Motion – 15 Days Negative Pressure Wound Therapy system

REF 087.4024 Invia Motion – 60 Days Negative Pressure Wound Therapy system

These two Invia Motion NPWT pump versions have a runtime of 15 and 60 days. After one hour of use, the timer is activated and will continue even when the pump is turned off. During the first hour of use, the Invia Motion NPWT pump can be turned on and off without starting the timer. If the timer function is used to readout the therapy time, one hour needs to be added to the displayed runtime. Press [] for information about the used operation time.

After 15 or 60 days, a blinking symbol on the display alike [] will indicate the expiration of the Invia Motion NPWT pump.

Example for 15 days version (blinking)

 The Invia Motion NPWT pump runtime versions of 15 or 60 days provide an additional safety usage time of 72 hours. During this time the device can be operated as usual. An acoustic signal will occur after 72 hours and the negative pressure will not be maintained.

Setup carrying case

The Invia Motion NPWT pump is intended to be worn in a carrying case during operation. An Invia Motion carrying case is provided with the Invia Motion NPWT pump.

Attach shoulder/belt strap onto the carrying case as a first step.

1. Open the carrying case.
2. Slide the Invia Motion NPWT pump into the designated upper pocket and close the Velcro strap.
3. Place the excess tubing into the designated lower pocket and close the Velcro strap.
4. Fold the pump and tubing sections of the carrying case together and align the inner Velcro tabs to secure.
5. Secure the bag by closing the lid and fastening the button.

The carrying case strap can be worn over the shoulder or around the waist.

To wear over the shoulder, attach the strap clips to the carrying case. To wear as a belt, thread the strap through the two loops located on the back of the case.

Shoulder strap

Belt strap

Operating instructions

Invia Motion NPWT pump settings

Change pressure level

Press [] and hold for 3 seconds to change the pressure level. The set pressure value symbol on the display will start to blink. You now have 5 seconds to press [] for a new pressure setting, otherwise the pump remains at the previous setting. Each time the button is pressed the next available pressure will appear. When the appropriate pressure level appears on display, the pump will adopt it automatically after 5 seconds. The pressure may only be changed if instructed by a physician or qualified healthcare professional.

| Available pressure levels | | | | | | | |
|---------------------------|-----|-----|-----|------|------|------|------|
| mmHg | -40 | -60 | -80 | -100 | -125 | -150 | -175 |

Change therapy mode


The Invia Motion NPWT pump is able to deliver constant or intermittent (5 min. on/2 min. off) suction.

Constant mode

Intermittent mode


Press [] and hold for 3 seconds to change the therapy mode. The therapy mode symbol on the display will start to blink. You now have 5 seconds to press [] to change the mode. Otherwise the pump remains at the previous mode. Each time the button is pressed the next available therapy mode will appear. When the appropriate mode symbol appears, the pump will adopt it automatically after 5 seconds. The therapy mode may only be changed if instructed by a physician or a qualified healthcare professional.

Pump in constant mode

Pump in intermittent mode (on)

Pump in intermittent mode (off)

Mute / Information

Press [] to obtain information about the runtime used or to mute the acoustic signal.

Change Invia Motion canister/tubing set

The Invia Motion canister/tubing set has to be changed on the basis of a visual check or according to the instructions on the display but at a minimum of once a week.

1. Wash hands and prepare clean paper towel.
2. Remove the Invia Motion NPWT pump from the carrying case. Switch the Invia Motion NPWT pump off by pressing [] and holding for 3 seconds. Press only once.
3. Position the clamp next to the quick-connector and close the clamp on the pump tubing and close the clamp on the dressing tubing.


Clamp open

Clamp closed

4. Disconnect the pump tubing from the dressing tubing by pressing on the sides of the quick-connector. Have a paper towel or tissue available to cover the tip of the used tubing to protect against wound fluid leakage.

5. Detach the full canister/tubing set by rotating the canister to the left.

6. Wash hands.
7. Open the sterile packaging of a new Invia Motion canister/tubing set.

CAUTION

Do not use an Invia Motion canister/tubing set if the sterile packaging is damaged or has expired.

Make sure the correct canister is used: **REF** 087.4016.

8. Remove the canister/tubing set from the packaging and attach it to Invia Motion NPWT pump.

For dressing applications, please refer to instructions for use provided with Invia dressing kits.

9. Connect the pump tubing to the dressing tubing. Push together until you hear a click.

 Make sure both parts are aligned correctly as shown in the picture below.

10. Open the clamp on the dressing tubing. Also make sure the clamp on the pump tubing is open.

11. Switch the Invia Motion NPWT pump on by pressing [] and slide it back into the carrying case.

12. Discard the used canister/tubing set according to local regulations.

Power supply and battery management

The Invia Motion NPWT pump is operated by an internal rechargeable NiMH battery and while attached to a power source. Make sure the battery compartment on the back side of the pump is closed. Do not remove battery from device as long as therapy is ongoing. If the rechargeable battery is removed more than 30 minutes the Invia Motion NPWT pump will suspend to operate.

For disposal of rechargeable battery after therapy, refer to chapter Disposal. Fully charge battery before first use. Charging time can be up to 4 hours.

While in use and connected to a power source, the battery is recharged. The available battery time from a single charge depends upon the runtime of the motor. This is influenced by the size of the wound, air leaks in the dressing, the set pressure and the mode of therapy. The pump only runs if the pressure drops below the set pressure. If the Invia Motion NPWT pump is running continuously, a maximum of 4 hours of battery time is expected before it is necessary to be recharged. The expected battery time generally exceeds 10 hours. If the Invia Motion NPWT pump is not in use, the battery must be charged approximately once every 6 months to ensure optimal functioning.

 Before attaching the Invia Motion NPWT pump to a power source, check to ensure your local power supply is the same as the voltage indicated on the specification plate located on the back of the pump.


Battery low

Approximately 30 minutes before the battery is fully discharged, an acoustic signal will sound and the battery symbol [] will start to blink. The acoustic signal can be muted, but the battery symbol continues to blink until the battery is fully discharged. The set pressure is maintained during the 30 minutes, but the battery should be recharged as soon as possible. Fully charge battery. Charging time can be up to 4 hours.

Battery fully discharged

The battery symbol [] blinks on the display, an acoustic signal sounds and the negative pressure is not maintained. Connect the pump to the power supply. When the pump is not connected to the power supply within 10 minutes, it switches off.

Connecting power source

Plug the Invia Motion power supply into a wall outlet and open the small protective cover to access the DC socket. Plug the coaxial power supply into the DC socket.

Moving bars on the battery display will indicate charging is occurring [].

1. Battery is fully discharged

2. Battery is being charged (moving bars)

3. Battery is fully charged

Disconnecting power source

When the battery is fully charged (after approximately 4 hours) the following non-blinking symbol will appear on the display []. Charging of the battery can be slowed when exposed to high temperatures. To disconnect the Invia Motion NPWT pump from the electrical source, hold the end of the plug and remove it from the pump.

Close the DC socket by pressing the small protective cover over the opening.

Notifications and troubleshooting

WARNING

Negative Pressure Wound Therapy should run 24 hours a day without interruption. If therapy is discontinued for more than 2 hours, the dressing should be replaced and therapy restarted by a healthcare professional.

If the Invia Motion NPWT pump detects an irregularity in the operation of the pump, an acoustic signal sounds and a notification will appear on the display.

1. An acoustic signal sounds.
2. If the Invia Motion NPWT pump is in the carrying case, open the lid of the carrying case and release the Velcro.


3. Press [] to mute the acoustic signal for 60 seconds (the button can be pressed every time the signal sounds, with the exception of "Battery fully discharged" and after the pump time of 15 or 60 days plus 72 hours of safety time is used up).

Every acoustic signal can be permanently muted by switching Invia Motion NPWT pump off (press [] and hold for 3 seconds). Turn the pump back on to continue treatment.

4. Check the display for one of the following visual symbols.

Canister full


Acoustic signal will sound when the Invia Motion NPWT pump detects a full canister and canister full symbol will appear on the display.

Troubleshooting:

Change the Invia Motion canister/tubing set.

Blockage in tubing

Acoustic signal will sound when the Invia Motion NPWT pump detects blockage in tubing and blockage symbol will appear on the display.

Troubleshooting:

Inspect Invia Motion canister and dressing tubing to ensure tubing is not twisted, kinked or clamped. Verify that the correct canister is used: **REF 087.4016**.

If blockage is identified along canister tubing, switch Invia Motion NPWT pump off and follow instructions in chapter "Operating instructions" to change the Invia Motion canister/tubing set. If the used canister is not the correct one, switch the Invia Motion NPWT pump off, remove the canister. Attach a new canister:

REF 087.4016.

If blockage is identified along dressing tubing, switch Invia Motion NPWT pump off and contact your healthcare professional to replace dressing and restart therapy.

If Invia Motion NPWT pump detects blockage in the intermittent mode, it will continue to run in the on-state until blockage is resolved.

The notification symbol and the acoustic signal may still be active, even if situation is resolved.

Wait a short time for the pump to react.


WARNING

Invia Motion pump does not detect blockage along tubing connected to the drain adapter or the curved arm of the Y-connector.

1. Instruct the patient to frequently check dressing tubing connected to the drain adapter or the curved arm of the Y-connector to make sure tubing is not kinked or blocked.
2. If blockage is identified along tubing connected to the drain adapter or tubing connected to the curved arm of the Y-connector, instruct the patient to contact his healthcare professional.


Leak in system

Acoustic signal will sound when the Invia Motion NPWT pump detects a leak and leakage symbol will appear on the display.

Troubleshooting:

Inspect dressing site for leaks and make sure the Invia Motion canister/tubing set quick-connector is fully engaged. Check that canister is secured to the pump. If a leak is identified at the dressing site, press down on dressing to seal.


The notification symbol and the acoustic signal may still be active, even if situation is resolved.

Wait for some time to let the pump react. If acoustic signal and leakage symbol do not disappear, switch Invia Motion NPWT Pump off and follow instructions in chapter "Operating instructions" to change the Invia Motion canister/tubing set.

☞ If leak is not resolved, switch Invia Motion NPWT pump off and contact your healthcare professional to replace dressing and restart therapy.

Battery low


Troubleshooting:

Charge the battery.

Approximately 30 minutes before the battery is fully discharged, an acoustic signal will sound and the battery symbol [] will start to blink. The acoustic signal can be muted, but the

battery symbol continues to blink until the battery is fully discharged. The set pressure is maintained, but the battery should be recharged as soon as possible. Fully charge battery. Charging time can be up to 4 hours.


Battery fully discharged


Troubleshooting:

Charge battery


The battery symbol [] blinks on the display, an acoustic signal sounds and the negative pressure is not maintained. Connect the pump to the power supply. When the pump is not connected to the power supply within 10 minutes, it switches off.


Info used time


The Invia Motion NPWT pumps have a runtime of 15 or 60 days. After one hour of use, the timer is activated and will continue even when the pump is turned off. During the first hour of use, the Invia Motion NPWT pump can be turned on and off without starting the timer. If the timer function is used to readout the therapy time, one hour needs to be added to the displayed runtime. Press [] for information about the used operation time.

The example below displays 6 days and 23 hours of operation time used. The remaining runtime of the Invia Motion – 15 Days version would therefore be 8 days and 1 hour.

Runtime expired



After 15 or 60 days, a blinking symbol on the display alike [] will indicate the expiration of the Invia Motion NPWT pump.

Example for 15 days version (blinking)

 The Invia Motion NPWT pump runtime versions of 15 or 60 days provide an additional safety usage time of 72 hours. During this time the device can be operated as usual. An acoustic signal will occur after 72 hours and the negative pressure will not be maintained.

Overview notifications

| | Visual symbol on display | Notification Description | Acoustic signal | Troubleshooting | Remarks/Potential cause |
|--------------|---|--------------------------|---|---|--|
| Notification |  Blinking! + acoustic signal! | Leak in system | Can be muted for 60 seconds by pressing  i | Check the system for leaks. | Dressing could be leaking, check that it is properly sealed. |
| |  Blinking! + acoustic signal! | Blockage | Can be muted for 60 seconds by pressing  i | Check the system for blockage. Verify that the correct canister is used: REF 087.4016. | The dressing or canister tubing could be twisted, kinked or clamped. If the Invia Motion pump detects blockage in the intermittent mode, it will continue to run in the on-state until blockage is resolved. |
| |  Blinking! + acoustic signal! | Canister full | Can be muted for 60 seconds by pressing  i | Change the Invia Motion canister/ tubing set. | The filter could get wet by improper positioning of the pump (upside down or display side down). |
| |  Blinking! + acoustic signal! | Battery low | Can be muted permanently by pressing  i | Connect the Invia Motion pump to the power supply and fully charge battery. Charging time can be up to 4 hours. | The visual symbol blinks all the time until battery is fully discharged. The remaining time of the battery is approximately 30 minutes. |
| |  Blinking! + acoustic signal! | Battery fully discharged | Can be muted by switching off the pump. | Connect the Invia Motion pump to the power supply and fully charge battery. Charging time can be up to 4 hours | The remaining time of the battery is approximately 10 minutes. |
| |  Blinking! (Example for 15 days version) | Pump time expired | The blinking signal can be suppressed permanently by pressing  i | Dispose of the Invia Motion NPWT pump. Attach another device if necessary. | The Invia Motion NPWT pump provides an additional safety usage time of 72 hours. During this time the Invia Motion pump can be operated as usual. A signal will occur after 72 hours and the negative pressure will not be maintained. |

Overview notifications

| | | | | | |
|-----------------------|--|-----------------|---|--|--|
| Notification |  Blinking! + acoustic signal! | Battery missing | Can be muted for 60 seconds by pressing  ; | Ensure that battery compartment cover is mounted and closed properly. | Do not remove battery during ongoing therapy. |
| Internal fault |  Acoustic signal! | Internal fault | Can be muted by switching off the pump. | Switch the Invia Motion NPWT pump off, release and refix canister/tubing set. Turn the Invia Motion pump back on. If the fault message repeats, contact Medela customer service. | If the fault message repeats, contact Medela customer service. |

Contact your local Medela customer service representative for assistance with product operations.

Cleaning and disinfection

Cleaning

WARNING

Unplug the pump before cleaning. Close the small protective cover on the DC socket.

The Invia Motion NPWT pump and carrying case are for a single patient's entire therapy. Periodic cleaning is recommended.

Wipe off with a damp cloth or use a non-abrasive detergent.

Disinfection

Invia Motion can be disinfected with the disinfecting agent group "alcohol".

Sterilization

The Invia Motion NPWT pump and accessories cannot be sterilized. The canister/tubing set is a disposable item and may not be reused or resterilized. Replace after every use!

The Invia Motion canister/tubing set is a single patient use product and is not intended to be reprocessed. Reprocessing may cause the loss of mechanical, chemical and/or biological characteristics.

WARNING

Reuse on another patient might cause cross contamination.

General notes

Comply with facility or institutional guidelines on hygiene and infection control. Dispose of fluids such as blood and fluids, along with all parts contaminated with them in accordance with facility requirements. Wear suitable gloves for cleaning/disinfection and disposal.

Disposal

Turn the Invia Motion NPWT pump off by pressing [] and holding for 3 seconds.

Disconnect the canister/tubing set from the pump.
Dispose of the canister/tubing set and pump separately and according to local regulations.

When the Invia Motion NPWT pump is ready for disposal, take care to dispose of the Invia Motion NPWT pump and its accessories in accordance with your local regulations.

Disposal of electronic equipment

The symbol above means that the electrical and electronic equipment must not be disposed as normal household refuse. The correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about proper disposal, contact Medela or your local healthcare professional. This symbol is only valid in the European Union. Please dispose according to specified institution protocol or state and local regulations of electrical and electronic equipment.

 Do not remove battery from device as long as therapy is ongoing.

Battery can be disposed separately after completion of therapy:

1. Cut the adhesive label.
2. Open battery cover.
3. Remove rechargeable battery.
4. Dispose rechargeable battery according to local regulations.

Warranty and maintenance

Warranty

The warranty period for the Invia Motion NPWT system expires after 15 or 60 days of usage plus 72 hours of safety time, if used in accordance with the instructions for use. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage or unauthorized persons using the device.

Maintenance

The Invia Motion NPWT pump is maintenance free and will not require service.

If an Invia Motion NPWT pump fails within the warranty period due to a manufacturing defect, the pump will be replaced. The original Invia Motion NPWT pump will need to be returned to the supplier.

Safety-related checks

The Invia Motion NPWT system is a device in electrical protection Class II. The safety-related checks are confined to visual inspection of the pump housing and power supply for damage. This check must be performed prior to each use.

Devices of protection Class II do not have a protective earth conductor. There is no need to check the earth leakage current.

The Invia Motion NPWT pump enclosure is made entirely of insulated material. Tests of the enclosure leakage current using common measuring instruments will therefore not reveal measurable values.

The Invia Motion NPWT pump does not have patient circuits or functional earth connections.


Signs and symbols

 This symbol indicates safety related tip.

 This symbol indicates a CAUTION or WARNING associated with the device (see page 7).

 This symbol indicates a class II device.

 This symbol indicates to not dispose the device together with household refuse (for EU only).

 This symbol indicates the date of manufacture.

 This symbol indicates the name and the address of the manufacturer.

 This symbol indicates the device is sterilized using ethylene oxide.

 This symbol indicates a prescription device. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician (for US only).

 This symbol indicates a type BF applied part.

 The symbol indicates the item is a medical device.

 This symbol indicates manufacturer's catalog number.

 This symbol indicates manufacturer's serial number.

 This symbol indicates manufacturer's batch code.

IP22 This symbol indicates the protection against the ingress of dripping water and solid foreign objects of 12,5mm and greater.

 This symbol indicates that the device should not be used after the date shown.

 This symbol indicates to follow the instruction for use.

 This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

 This symbol indicates to not use the device if package is damaged.

pcs This symbol indicates the number of items

 This symbol indicates the direct current 9 VDC socket.

This symbol indicates a single use device.
Do not reuse the device.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates MR unsafe.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates to keep the device dry.

This symbol indicates to handle the fragile device with care.

This symbol indicates to keep the device away from sunlight.

This symbol indicates that the marked item or its material is part of a recovery or recycling process.

This symbol indicates the correct upright position of the transport package.

This symbol indicates that the device is in conformance with the Medical Device Directive 93/42/EEC.

Technical specifications

medium vacuum
-175 mmHg

low flow

350 g
Without canister


IP22


H x W x D
99 x 98 x 52 mm

ISO 9001
ISO 13485
CE (93/42/EEC), IIa

Operating
conditions

Transport/
Storage
conditions

Invia Motion power supply

Model:

UES12LCP-090133SPA

Input: 100-240VAC,
0.5A, 50-60 Hz

Output: 9VDC, 1.33A

Length of power supply
cable: 180 cm

Invia Motion Pump

[VDC] 9.0

[W] 8.0

Invia Motion battery

Type: NiMH

Invia Motion canister/tubing set

Material: Polypropylene, Silicon, ABS

Sterilization Method: Ethylene Oxide

Electromagnetic compatibility (EMC)

Invia Motion is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according clause 7 and 8.9.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition and IEC 60601-1-2:2014 4th Edition, Table 1)

| Electromagnetic Emissions | | |
|---|-------------------|---|
| The Invia Motion NPWT system is intended for use in the electromagnetic environment specified below. The customer or the user of Invia Motion NPWT system should assure that they are used in such environment. | | |
| Emission Tests | Compliance | Electromagnetic environment – guidance |
| RF emissions CISPR 11 | Group 1 | The Invia Motion NPWT system uses RF energy only for their internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations / flicker emissions IEC 61000-3-3 | Complies | |

WARNING

The Invia Motion NPWT system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Motion NPWT system should be observed to verify normal operation in the configuration in which it will be used.

WARNING

Invia Motion is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

WARNING

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Invia Motion pump.

WARNING

EMC

HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Invia Motion pump.

WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Invia Motion. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic compatibility
(EMC, IEC 60601-1-2:2014 4th Edition, Table 2)

| Electromagnetic immunity | | | |
|--|--|--|--|
| The Invia Motion pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Motion pump should assure that it is used in such environment. | | | |
| Immunity tests | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
| Electrostatic Discharge (ESD) IEC 61000-4-2 | ± 8 kV contact ± 15 kV air | ± 8 kV contact ± 15 kV air | The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input / output lines | ± 2 kV for power supply lines ± 1 kV for input / output lines not applicable | The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Surge IEC 61000-4-5 | ± 1 kV differential mode ± 2 kV line-to-earth | ± 1 kV differential mode ± 2 kV line-to-earth not applicable | The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0% U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T for 1 cycle 70% U_T for 25 cycles at 50 Hz Single phase: at 0° 0% U_T for 250 cycles at 50 Hz | 0% U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T for 1 cycle 70% U_T for 25 cycles at 50 Hz Single phase: at 0° 0% U_T for 250 cycles at 50 Hz | The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. If the user of the Invia Motion pump requires continued operation during power mains interruptions, it is recommended that the Invia Motion pump be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m | 30 A/m | The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| NOTE: U_T is the a.c. mains voltage prior to application of the test level. | | | |

Electromagnetic compatibility
(EMC, IEC 60601-1-2:2014 4th Edition, Table 4)

| Electromagnetic immunity | | | |
|---|---|-------------------------|---|
| <p>The Invia Motion pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Motion pump should assure that it is used in such environment.</p> | | | |
| Immunity tests | IEC 60601 test level | Compliance level | Electromagnetic environment-guidance |
| <p>Conducted RF IEC 61000-4-6</p> | <p>3 Vrms 150 kHz to 80 MHz</p> | <p>3 Vrms</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the Invia Motion pump including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ $d = 0.35 \sqrt{P}$ <p style="text-align: right;">80 MHz to 800 MHz</p> $d = 0.7 \sqrt{P}$ <p style="text-align: right;">800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| <p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> | | | |
| <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> | | | |
| <p>^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Invia Motion pump is used exceeds the applicable RF compliance level above, the Invia Motion pump should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating the Invia Motion pump.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.</p> | | | |

Electromagnetic compatibility
 (EMC, IEC 60601-1-2:2014 4th Edition, Table 6)

Recommended separation distance between portable and mobile RF communications equipment and the Invia Motion pump

The Invia Motion pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Invia Motion pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Invia Motion pump as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter | Separation distance according to frequency of transmitter | | |
|---|---|---|---|
| | m | | |
| W | 150 kHz to 80 MHz $d = 1.2 \sqrt{P}$ | 80 MHz to 800 MHz $d = 0.35 \sqrt{P}$ | 800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$ |
| 0.01 | 0.12 | 0.035 | 0.07 |
| 0.1 | 0.38 | 0.11 | 0.22 |
| 1 | 1.2 | 0.35 | 0.7 |
| 10 | 3.8 | 1.11 | 2.21 |
| 100 | 12 | 3.5 | 7.0 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance (*d*) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Medela AG
Medela UK Ltd.